

ALK Conjugated Antibody

Catalog No: #C56402

Package Size: #C56402-AF350 100ul #C56402-AF405 100ul #C56402-AF488 100ul #C56402-AF555 100ul #C56402-AF594 100ul #C56402-AF647 100ul #C56402-AF680 100ul #C56402-AF750 100ul #C56402-Biotin 100ul #C56402-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ALK Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IF
Species Reactivity	Human Mouse Rat
Specificity	ALK Antibody detects endogenous levels of total ALK
Immunogen Description	A synthesized peptide derived from human ALK
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ALK tyrosine kinase receptor; Nucleophosmin; NPM1; ALK; NPM-ALK
Accession No.	Uniprot:Q9UM73
Calculated MW	kDa
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

Product Description

Involved in diverse cellular processes such as ribosome biogenesis, centrosome duplication, protein chaperoning, histone assembly, cell proliferation, and regulation of tumor suppressors p53/TP53 and ARF. Binds ribosome presumably to drive ribosome nuclear export. Associated with nucleolar ribonucleoprotein structures and bind single-stranded nucleic acids.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.