

PBRM1 Conjugated Antibody

Catalog No: #C56406

Package Size: #C56406-AF350 100ul #C56406-AF405 100ul #C56406-AF488 100ul #C56406-AF555 100ul #C56406-AF594 100ul #C56406-AF647 100ul #C56406-AF680 100ul #C56406-AF750 100ul #C56406-Biotin 100ul #C56406-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PBRM1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PBRM1 Antibody detects endogenous levels of total PBRM1
Immunogen Description	A synthesized peptide derived from human PBRM1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BAF180; hPB1; PB1; Pbrm1; polybromo 1; Polybromo-1D;
Accession No.	Uniprot:Q86U86
Calculated MW	193kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology). Act as a negative regulator of cell proliferation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.