

Parathyroid Hormone Receptor 1 Conjugated Antibody

Catalog No: #C56417

Orders: order@signalwayantibody.com

Package Size: #C56417-AF350 100ul #C56417-AF405 100ul #C56417-AF488 100ul #C56417-AF555 100ul #C56417-AF594 100ul

#C56417-AF647 100ul #C56417-AF680 100ul #C56417-AF750 100ul #C56417-Biotin 100ul #C56417-Conjugated 50ul

Description

Product Name	Parathyroid Hormone Receptor 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Parathyroid Hormone Receptor 1 Antibody detects endogenous levels of total Parathyroid Hormone Receptor 1
Immunogen Description	A synthesized peptide derived from human Parathyroid Hormone Receptor 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PTH receptor; PTH1 receptor; PTH1R; PTHR 1; PTHR; PTHR1;
Accession No.	Uniprot:Q03431
Calculated MW	66kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Receptor for parathyroid hormone and for parathyroid hormone-related peptide. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase and also a phosphatidylinositol-calcium second messenger system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.