

# CD134 Conjugated Antibody

Catalog No: #C56422

Package Size: #C56422-AF350 100ul #C56422-AF405 100ul #C56422-AF488 100ul #C56422-AF555 100ul #C56422-AF594 100ul #C56422-AF647 100ul #C56422-AF680 100ul #C56422-AF750 100ul #C56422-Biotin 100ul #C56422-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD134 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	CD134 Antibody detects endogenous levels of total CD134
Immunogen Description	A synthesized peptide derived from human CD134
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ACT35; CD134; IMD16; OX40; OX40L receptor; TNFRSF4; Txgp1; Txgp1l;
Accession No.	Uniprot:P43489
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Receptor for TNFSF4/OX40/GP34.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.