

GSTO1 Conjugated Antibody

Catalog No: #C56435

Package Size: #C56435-AF350 100ul #C56435-AF405 100ul #C56435-AF488 100ul #C56435-AF555 100ul #C56435-AF594 100ul #C56435-AF647 100ul #C56435-AF680 100ul #C56435-AF750 100ul #C56435-Biotin 100ul #C56435-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GSTO1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Rat
Specificity	GSTO1 Antibody detects endogenous levels of total GSTO1
Immunogen Description	A synthesized peptide derived from human GSTO1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GSTO1 1; Gsto1; GSTTLp28; GSTX; Gtsttl; p28; SPG-R;
Accession No.	Uniprot:P78417
Calculated MW	28kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Exhibits glutathione-dependent thiol transferase and dehydroascorbate reductase activities. Has S-(phenacyl)glutathione reductase activity. Has also glutathione S-transferase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.