

PSMD4 Conjugated Antibody

Catalog No: #C56439

Package Size: #C56439-AF350 100ul #C56439-AF405 100ul #C56439-AF488 100ul #C56439-AF555 100ul #C56439-AF594 100ul #C56439-AF647 100ul #C56439-AF680 100ul #C56439-AF750 100ul #C56439-Biotin 100ul #C56439-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PSMD4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PSMD4 Antibody detects endogenous levels of total PSMD4
Immunogen Description	A synthesized peptide derived from human PSMD4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PSMD4; AF; AF-1; ASF; MCB1; Rpn10; S5A; pUB-R5;
Accession No.	Uniprot:P55036
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Binds and presumably selects ubiquitin-conjugates for destruction. Displays selectivity for longer polyubiquitin chains. Modulates intestinal fluid secretion.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.