

eIF4A1 Conjugated Antibody

Catalog No: #C56446

Package Size: #C56446-AF350 100ul #C56446-AF405 100ul #C56446-AF488 100ul #C56446-AF555 100ul #C56446-AF594 100ul #C56446-AF647 100ul #C56446-AF680 100ul #C56446-AF750 100ul #C56446-Biotin 100ul #C56446-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	eIF4A1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	eIF4A1 Antibody detects endogenous levels of total eIF4A1
Immunogen Description	A synthesized peptide derived from human eIF4A1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DDX2; DDX2A; eIF 4A I; EIF4A; eIF4A I; eIF4A-I; EIF4A1;
Accession No.	Uniprot:P60842
Calculated MW	46kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

ATP-dependent RNA helicase which is a subunit of the eIF4F complex involved in cap recognition and is required for mRNA binding to ribosome.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.