

IRS2 Conjugated Antibody

Catalog No: #C56455

Package Size: #C56455-AF350 100ul #C56455-AF405 100ul #C56455-AF488 100ul #C56455-AF555 100ul #C56455-AF594 100ul #C56455-AF647 100ul #C56455-AF680 100ul #C56455-AF750 100ul #C56455-Biotin 100ul #C56455-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IRS2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	IRS2 Antibody detects endogenous levels of total IRS2
Immunogen Description	A synthesized peptide derived from human IRS2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Insulin receptor substrate 2; IRS 2; IRS-2;
Accession No.	Uniprot:Q9Y4H2
Calculated MW	185kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

May mediate the control of various cellular processes by insulin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.