

Prostatic Acid Phosphatase Conjugated Antibody

Catalog No: #C56459



Package Size: #C56459-AF350 100ul #C56459-AF405 100ul #C56459-AF488 100ul #C56459-AF555 100ul #C56459-AF594 100ul #C56459-AF647 100ul #C56459-AF680 100ul #C56459-AF750 100ul #C56459-Biotin 100ul #C56459-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Prostatic Acid Phosphatase Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse
Specificity	Prostatic Acid Phosphatase Antibody detects endogenous levels of total Prostatic Acid Phosphatase
Immunogen Description	A synthesized peptide derived from human Prostatic Acid Phosphatase
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ACP3; PAP; TMPase;
Accession No.	Uniprot:P15309
Calculated MW	52kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

A non-specific tyrosine phosphatase that dephosphorylates a diverse number of substrates under acidic conditions (pH 4-6) including alkyl, aryl, and acyl orthophosphate monoesters and phosphorylated proteins. Has lipid phosphatase activity and inactivates lysophosphatidic acid in seminal plasma.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.