

PIAS1 Conjugated Antibody

Catalog No: #C56461

Package Size: #C56461-AF350 100ul #C56461-AF405 100ul #C56461-AF488 100ul #C56461-AF555 100ul #C56461-AF594 100ul #C56461-AF647 100ul #C56461-AF680 100ul #C56461-AF750 100ul #C56461-Biotin 100ul #C56461-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PIAS1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PIAS1 Antibody detects endogenous levels of total PIAS1
Immunogen Description	A synthesized peptide derived from human PIAS1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DDXBP1; GBP; GU/RH-II; Pias1; ZMIZ3;
Accession No.	Uniprot:O75925
Calculated MW	76kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Functions as an E3-type small ubiquitin-like modifier (SUMO) ligase, stabilizing the interaction between UBE2I and the substrate, and as a SUMO-tethering factor. Plays a crucial role as a transcriptional coregulation in various cellular pathways, including the STAT pathway, the p53 pathway and the steroid hormone signaling pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.