

Eph receptor B1 Conjugated Antibody

Catalog No: #C56463

Package Size: #C56463-AF350 100ul #C56463-AF405 100ul #C56463-AF488 100ul #C56463-AF555 100ul #C56463-AF594 100ul #C56463-AF647 100ul #C56463-AF680 100ul #C56463-AF750 100ul #C56463-Biotin 100ul #C56463-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Eph receptor B1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Eph receptor B1 Antibody detects endogenous levels of total Eph receptor B1
Immunogen Description	A synthesized peptide derived from human Eph receptor B1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ELK; NET; Hek6; EPHT2; EPHB1;
Accession No.	Uniprot:P54762
Calculated MW	120kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Receptor for members of the ephrin-B family. Binds to ephrin-B1, -B2 and -B3. Binding with the guidance cue ephrin-B2 at the optic chiasm midline redirect ventrotemporal (VT) retinal ganglion cells (RGCs) axons ipsilaterally. May be involved in cell-cell interactions in the nervous system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.