

# Cripto1 Conjugated Antibody

Catalog No: #C56545

Package Size: #C56545-AF350 100ul #C56545-AF405 100ul #C56545-AF488 100ul #C56545-AF555 100ul #C56545-AF594 100ul #C56545-AF647 100ul #C56545-AF680 100ul #C56545-AF750 100ul #C56545-Biotin 100ul #C56545-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cripto1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Cripto1 Antibody detects endogenous levels of total Cripto1
Immunogen Description	A synthesized peptide derived from human Cripto1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CRGF; Cripto 1; cripto; Cripto1 growth factor; TDGF1;
Accession No.	Uniprot:P13385
Calculated MW	21 kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Could play a role in the determination of the epiblastic cells that subsequently give rise to the mesoderm.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.