

## GPA33 Conjugated Antibody

Catalog No: #C56551



Package Size: #C56551-AF350 100ul #C56551-AF405 100ul #C56551-AF488 100ul #C56551-AF555 100ul #C56551-AF594 100ul

#C56551-AF647 100ul #C56551-AF680 100ul #C56551-AF750 100ul #C56551-Biotin 100ul #C56551-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPA33 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	GPA33 Antibody detects endogenous levels of total GPA33
Immunogen Description	A synthesized peptide derived from human GPA33
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Cell surface A33 antigen; Gpa33;
Accession No.	Uniprot:Q99795
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

May play a role in cell-cell recognition and signaling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.