

# CD16 Conjugated Antibody

Catalog No: #C56554

Package Size: #C56554-AF350 100ul #C56554-AF405 100ul #C56554-AF488 100ul #C56554-AF555 100ul #C56554-AF594 100ul #C56554-AF647 100ul #C56554-AF680 100ul #C56554-AF750 100ul #C56554-Biotin 100ul #C56554-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD16 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CD16 Antibody detects endogenous levels of total CD16
Immunogen Description	A synthesized peptide derived from human CD16
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	IGFR3; CD16; CD16a; IMD20;
Accession No.	Uniprot:P08637
Calculated MW	42kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Receptor for the Fc region of IgG. Binds complexed or aggregated IgG and also monomeric IgG. Mediates antibody-dependent cellular cytotoxicity (ADCC) and other antibody-dependent responses, such as phagocytosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.