

STK33 Conjugated Antibody

Catalog No: #C56567

Package Size: #C56567-AF350 100ul #C56567-AF405 100ul #C56567-AF488 100ul #C56567-AF555 100ul #C56567-AF594 100ul #C56567-AF647 100ul #C56567-AF680 100ul #C56567-AF750 100ul #C56567-Biotin 100ul #C56567-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	STK33 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	STK33 Antibody detects endogenous levels of total STK33
Immunogen Description	A synthesized peptide derived from human STK33
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Stk33;
Accession No.	Uniprot:Q9BYT3
Calculated MW	58, 50kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Serine/threonine protein kinase which phosphorylates VIME. May play a specific role in the dynamic behavior of the intermediate filament cytoskeleton by phosphorylation of VIME (By similarity). Not essential for the survival of KRAS-dependent AML cell lines.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.