

# CD97 Conjugated Antibody

Catalog No: #C56568

Package Size: #C56568-AF350 100ul #C56568-AF405 100ul #C56568-AF488 100ul #C56568-AF555 100ul #C56568-AF594 100ul #C56568-AF647 100ul #C56568-AF680 100ul #C56568-AF750 100ul #C56568-Biotin 100ul #C56568-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD97 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Rat
Specificity	CD97 Antibody detects endogenous levels of total CD97
Immunogen Description	A synthesized peptide derived from human CD97
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ADGRE5; CD97; TM7LN1;
Accession No.	Uniprot:P48960
Calculated MW	55/28kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.