

GSC Conjugated Antibody

Catalog No: #C56573



Package Size: #C56573-AF350 100ul #C56573-AF405 100ul #C56573-AF488 100ul #C56573-AF555 100ul #C56573-AF594 100ul #C56573-AF647 100ul #C56573-AF680 100ul #C56573-AF750 100ul #C56573-Biotin 100ul #C56573-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GSC Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	GSC Antibody detects endogenous levels of total GSC
Immunogen Description	A synthesized peptide derived from human GSC
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Goosecoid; GSC;
Accession No.	Uniprot:P56915
Calculated MW	28kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Regulates chordin (CHRD). May play a role in spatial programming within discrete embryonic fields or lineage compartments during organogenesis. In concert with NKX3-2, plays a role in defining the structural components of the middle ear; required for the development of the entire tympanic ring.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.