

# BMP15 Conjugated Antibody

Catalog No: #C56576

Package Size: #C56576-AF350 100ul #C56576-AF405 100ul #C56576-AF488 100ul #C56576-AF555 100ul #C56576-AF594 100ul #C56576-AF647 100ul #C56576-AF680 100ul #C56576-AF750 100ul #C56576-Biotin 100ul #C56576-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	BMP15 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	BMP15 Antibody detects endogenous levels of total BMP15
Immunogen Description	A synthesized peptide derived from human BMP15
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BMP15; GDF9B; Growth/differentiation factor 9B; ODG2; POF4;
Accession No.	Uniprot:O95972
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

May be involved in follicular development. Oocyte-specific growth/differentiation factor that stimulates folliculogenesis and granulosa cell (GC) growth.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.