

PSMA2 Conjugated Antibody

Catalog No: #C56577

Package Size: #C56577-AF350 100ul #C56577-AF405 100ul #C56577-AF488 100ul #C56577-AF555 100ul #C56577-AF594 100ul #C56577-AF647 100ul #C56577-AF680 100ul #C56577-AF750 100ul #C56577-Biotin 100ul #C56577-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PSMA2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PSMA2 Antibody detects endogenous levels of total PSMA2
Immunogen Description	A synthesized peptide derived from human PSMA2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HC3; PMSA2; Proteasome component C3; PSC2; PSC3; psmA2;
Accession No.	Uniprot:P25787
Calculated MW	25kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.