

TIMM50 Conjugated Antibody

Catalog No: #C56579

Package Size: #C56579-AF350 100ul #C56579-AF405 100ul #C56579-AF488 100ul #C56579-AF555 100ul #C56579-AF594 100ul #C56579-AF647 100ul #C56579-AF680 100ul #C56579-AF750 100ul #C56579-Biotin 100ul #C56579-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TIMM50 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse
Specificity	TIMM50 Antibody detects endogenous levels of total TIMM50
Immunogen Description	A synthesized peptide derived from human TIMM50
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TIM50; TIM50L; TIMM50;
Accession No.	Uniprot:Q3ZCQ8
Calculated MW	36kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Essential component of the TIM23 complex, a complex that mediates the translocation of transit peptide-containing proteins across the mitochondrial inner membrane. Has some phosphatase activity in vitro; however such activity may not be relevant in vivo.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.