

## BCL2L12 Conjugated Antibody

Catalog No: #C56580



Package Size: #C56580-AF350 100ul #C56580-AF405 100ul #C56580-AF488 100ul #C56580-AF555 100ul #C56580-AF594 100ul #C56580-AF647 100ul #C56580-AF680 100ul #C56580-AF750 100ul #C56580-Biotin 100ul #C56580-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	BCL2L12 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	BCL2L12 Antibody detects endogenous levels of total BCL2L12
Immunogen Description	A synthesized peptide derived from human BCL2L12
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BCL2 like 12; BCL2L12; BPR;
Accession No.	Uniprot:Q9HB09
Calculated MW	32kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Belongs to the Bcl-2 family.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.