

WAPL Conjugated Antibody

Catalog No: #C56582

Package Size: #C56582-AF350 100ul #C56582-AF405 100ul #C56582-AF488 100ul #C56582-AF555 100ul #C56582-AF594 100ul #C56582-AF647 100ul #C56582-AF680 100ul #C56582-AF750 100ul #C56582-Biotin 100ul #C56582-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	WAPL Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	WAPL Antibody detects endogenous levels of total WAPL
Immunogen Description	A synthesized peptide derived from human WAPL
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	FOE; WAPAL; WAPL; Wings apart like homolog;
Accession No.	Uniprot:Q7Z5K2
Calculated MW	160kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Regulator of sister chromatid cohesion in mitosis which negatively regulates cohesin association with chromatin. Involved in both sister chromatid cohesion during interphase and sister-chromatid resolution during early stages of mitosis. Couples DNA replication to sister chromatid cohesion.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.