

Inhibin alpha Conjugated Antibody

Catalog No: #C56586

Package Size: #C56586-AF350 100ul #C56586-AF405 100ul #C56586-AF488 100ul #C56586-AF555 100ul #C56586-AF594 100ul #C56586-AF647 100ul #C56586-AF680 100ul #C56586-AF750 100ul #C56586-Biotin 100ul #C56586-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Inhibin alpha Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Inhibin alpha Antibody detects endogenous levels of total Inhibin alpha
Immunogen Description	A synthesized peptide derived from human Inhibin alpha
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	IHA; inhA; Inhibin alpha chain;
Accession No.	Uniprot:P05111
Calculated MW	40kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.