

SPINK1 Conjugated Antibody

Catalog No: #C56587



Package Size: #C56587-AF350 100ul #C56587-AF405 100ul #C56587-AF488 100ul #C56587-AF555 100ul #C56587-AF594 100ul #C56587-AF647 100ul #C56587-AF680 100ul #C56587-AF750 100ul #C56587-Biotin 100ul #C56587-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	SPINK1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	SPINK1 Antibody detects endogenous levels of total SPINK1
Immunogen Description	A synthesized peptide derived from human SPINK1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PCTT; PSTI; SPINK1; Spink3; TaTI; TCP;
Accession No.	Uniprot:P00995
Calculated MW	11kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Serine protease inhibitor which exhibits anti-trypsin activity (PubMed:7142173). In the pancreas, protects against trypsin-catalyzed premature activation of zymogens.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.