

## Nogo Receptor Conjugated Antibody

Catalog No: #C56600



Package Size: #C56600-AF350 100ul #C56600-AF405 100ul #C56600-AF488 100ul #C56600-AF555 100ul #C56600-AF594 100ul #C56600-AF647 100ul #C56600-AF680 100ul #C56600-AF750 100ul #C56600-Biotin 100ul #C56600-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nogo Receptor Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse
Specificity	Nogo Receptor Antibody detects endogenous levels of total Nogo Receptor
Immunogen Description	A synthesized peptide derived from human Nogo Receptor
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	NgR; Nogo 66 receptor; Nogo receptor; Nogor; Reticulon 4 receptor; Rtn4r;
Accession No.	Uniprot:Q9BZR6
Calculated MW	66kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Nogo Receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system. It is also a receptor for Nogo A (RTN4), OMG and MAG.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.