

Casein Conjugated Antibody

Catalog No: #C56603

Package Size: #C56603-AF350 100ul #C56603-AF405 100ul #C56603-AF488 100ul #C56603-AF555 100ul #C56603-AF594 100ul #C56603-AF647 100ul #C56603-AF680 100ul #C56603-AF750 100ul #C56603-Biotin 100ul #C56603-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Casein Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Casein Antibody detects endogenous levels of total Casein
Immunogen Description	A synthesized peptide derived from human Casein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Casein;
Accession No.	Uniprot:P47710
Calculated MW	22kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Important role in the capacity of milk to transport calcium phosphate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.