

TCF21 Conjugated Antibody

Catalog No: #C56612

Package Size: #C56612-AF350 100ul #C56612-AF405 100ul #C56612-AF488 100ul #C56612-AF555 100ul #C56612-AF594 100ul #C56612-AF647 100ul #C56612-AF680 100ul #C56612-AF750 100ul #C56612-Biotin 100ul #C56612-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TCF21 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	TCF21 Antibody detects endogenous levels of total TCF21
Immunogen Description	A synthesized peptide derived from human TCF21
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Capsulin; Epicardin; POD 1; POD1; TCF-21; tcf21;
Accession No.	Uniprot:O43680
Calculated MW	26kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Involved in epithelial-mesenchymal interactions in kidney and lung morphogenesis that include epithelial differentiation and branching morphogenesis. May play a role in the specification or differentiation of one or more subsets of epicardial cell types.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.