

Bmf Conjugated Antibody

Catalog No: #C56614

Package Size: #C56614-AF350 100ul #C56614-AF405 100ul #C56614-AF488 100ul #C56614-AF555 100ul #C56614-AF594 100ul #C56614-AF647 100ul #C56614-AF680 100ul #C56614-AF750 100ul #C56614-Biotin 100ul #C56614-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Bmf Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Bmf Antibody detects endogenous levels of total Bmf
Immunogen Description	A synthesized peptide derived from human Bmf
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Bcl2 modifying factor; Bmf;
Accession No.	Uniprot:Q96LC9
Calculated MW	24kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

May play a role in apoptosis. Isoform 1 seems to be the main initiator.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.