

VPS11 Conjugated Antibody

Catalog No: #C56615

Package Size: #C56615-AF350 100ul #C56615-AF405 100ul #C56615-AF488 100ul #C56615-AF555 100ul #C56615-AF594 100ul #C56615-AF647 100ul #C56615-AF680 100ul #C56615-AF750 100ul #C56615-Biotin 100ul #C56615-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	VPS11 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	VPS11 Antibody detects endogenous levels of total VPS11
Immunogen Description	A synthesized peptide derived from human VPS11
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	END1; hVPS11; PEP5; PP3476; RNF108; vps11;
Accession No.	Uniprot:Q9H270
Calculated MW	108kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

May play a role in vesicle-mediated protein trafficking to lysosomal compartments and in membrane docking/fusion reactions of late endosomes/lysosomes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.