

Y14 Conjugated Antibody

Catalog No: #C56616

Package Size: #C56616-AF350 100ul #C56616-AF405 100ul #C56616-AF488 100ul #C56616-AF555 100ul #C56616-AF594 100ul #C56616-AF647 100ul #C56616-AF680 100ul #C56616-AF750 100ul #C56616-Biotin 100ul #C56616-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Y14 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Y14 Antibody detects endogenous levels of total Y14
Immunogen Description	A synthesized peptide derived from human Y14
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BOV1; BOV1A; BOV1B; BOV1C; HSPC114; MDS014; RBM 8; RBM 8A; RBM 8B; RBM8; rbm8a; RBM8B; ZNRP; ZRNP1;
Accession No.	Uniprot:Q9Y5S9
Calculated MW	19kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

Component of a splicing-dependent multiprotein exon junction complex (EJC) deposited at splice junction on mRNAs. The EJC is a dynamic structure consisting of a few core proteins and several more peripheral nuclear and cytoplasmic associated factors that join the complex only transiently either during EJC assembly or during subsequent mRNA metabolism.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.