

ABCE1 Conjugated Antibody

Catalog No: #C56617



Package Size: #C56617-AF350 100ul #C56617-AF405 100ul #C56617-AF488 100ul #C56617-AF555 100ul #C56617-AF594 100ul #C56617-AF647 100ul #C56617-AF680 100ul #C56617-AF750 100ul #C56617-Biotin 100ul #C56617-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ABCE1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	ABCE1 Antibody detects endogenous levels of total ABCE1
Immunogen Description	A synthesized peptide derived from human ABCE1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ABC38; ABCE1; HuHP68; OABP; RLI; RNASEL1; RNASELI; RNS4I;
Accession No.	Uniprot:P61221
Calculated MW	67kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Antagonizes the binding of 2-5A (5'-phosphorylated 2',5'-linked oligoadenylates) by RNase L through direct interaction with RNase L and therefore inhibits its endoribonuclease activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.