

VPS18 Conjugated Antibody

Catalog No: #C56621

Package Size: #C56621-AF350 100ul #C56621-AF405 100ul #C56621-AF488 100ul #C56621-AF555 100ul #C56621-AF594 100ul #C56621-AF647 100ul #C56621-AF680 100ul #C56621-AF750 100ul #C56621-Biotin 100ul #C56621-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	VPS18 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	VPS18 Antibody detects endogenous levels of total VPS18
Immunogen Description	A synthesized peptide derived from human VPS18
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	hVPS18; PEP3; vps18;
Accession No.	Uniprot:Q9P253
Calculated MW	100kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

May play a role in vesicle-mediated protein trafficking to lysosomal compartments and in membrane docking/fusion reactions of late endosomes/lysosomes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.