

POLR2C Conjugated Antibody

Catalog No: #C56622

Package Size: #C56622-AF350 100ul #C56622-AF405 100ul #C56622-AF488 100ul #C56622-AF555 100ul #C56622-AF594 100ul #C56622-AF647 100ul #C56622-AF680 100ul #C56622-AF750 100ul #C56622-Biotin 100ul #C56622-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	POLR2C Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	POLR2C Antibody detects endogenous levels of total POLR2C
Immunogen Description	A synthesized peptide derived from human POLR2C
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	hRPB33; hsRPB3; Polr2c; RPB3; RPB31; RPB33;
Accession No.	Uniprot:P19387
Calculated MW	33kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Component of RNA polymerase II which synthesizes mRNA precursors and many functional non-coding RNAs.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.