

MYL9 Conjugated Antibody

Catalog No: #C56625

Package Size: #C56625-AF350 100ul #C56625-AF405 100ul #C56625-AF488 100ul #C56625-AF555 100ul #C56625-AF594 100ul #C56625-AF647 100ul #C56625-AF680 100ul #C56625-AF750 100ul #C56625-Biotin 100ul #C56625-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MYL9 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	MYL9 Antibody detects endogenous levels of total MYL9
Immunogen Description	A synthesized peptide derived from human MYL9
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MLC2; CMH10;
Accession No.	Uniprot:P24844
Calculated MW	18kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Myosin regulatory subunit that plays an important role in regulation of both smooth muscle and nonmuscle cell contractile activity via its phosphorylation. Implicated in cytokinesis, receptor capping, and cell locomotion.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.