

PAR2 Conjugated Antibody

Catalog No: #C56626

Package Size: #C56626-AF350 100ul #C56626-AF405 100ul #C56626-AF488 100ul #C56626-AF555 100ul #C56626-AF594 100ul #C56626-AF647 100ul #C56626-AF680 100ul #C56626-AF750 100ul #C56626-Biotin 100ul #C56626-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PAR2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PAR2 Antibody detects endogenous levels of total PAR2
Immunogen Description	A synthesized peptide derived from human PAR2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	F2RL1; GPR11; PAR 2;
Accession No.	Uniprot:P55085
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Receptor for trypsin and trypsin-like enzymes coupled to G proteins that stimulate phosphoinositide hydrolysis. May have a role in the regulation of vascular tone.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.