

# VAV3 Conjugated Antibody

Catalog No: #C56652

Package Size: #C56652-AF350 100ul #C56652-AF405 100ul #C56652-AF488 100ul #C56652-AF555 100ul #C56652-AF594 100ul #C56652-AF647 100ul #C56652-AF680 100ul #C56652-AF750 100ul #C56652-Biotin 100ul #C56652-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	VAV3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	VAV3 Antibody detects endogenous levels of total VAV3
Immunogen Description	A synthesized peptide derived from human VAV3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	RGD1565941; VAV 3; Vav3; VAV3 oncogene;
Accession No.	Uniprot:Q9UKW4
Calculated MW	98kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Plays an important role in angiogenesis. Its recruitment by phosphorylated EPHA2 is critical for EFNA1-induced RAC1 GTPase activation and vascular endothelial cell migration and assembly.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.