

# PDZK1 Conjugated Antibody

Catalog No: #C56686



Package Size: #C56686-AF350 100ul #C56686-AF405 100ul #C56686-AF488 100ul #C56686-AF555 100ul #C56686-AF594 100ul #C56686-AF647 100ul #C56686-AF680 100ul #C56686-AF750 100ul #C56686-Biotin 100ul #C56686-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PDZK1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	PDZK1 Antibody detects endogenous levels of total PDZK1
Immunogen Description	A synthesized peptide derived from human PDZK1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CAP70; CLAMP; PDZD1; NHERF3; NHERF-3;
Accession No.	Uniprot:Q5T2W1
Calculated MW	70kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

A scaffold protein that connects plasma membrane proteins and regulatory components, regulating their surface expression in epithelial cells apical domains.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.