

Pumilio 2 Conjugated Antibody

Catalog No: #C56688

Package Size: #C56688-AF350 100ul #C56688-AF405 100ul #C56688-AF488 100ul #C56688-AF555 100ul #C56688-AF594 100ul #C56688-AF647 100ul #C56688-AF680 100ul #C56688-AF750 100ul #C56688-Biotin 100ul #C56688-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Pumilio 2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Pumilio 2 Antibody detects endogenous levels of total Pumilio 2
Immunogen Description	A synthesized peptide derived from human Pumilio 2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PUM2; PUMH2; Pumilio2; PUML2;
Accession No.	Uniprot:Q8TB72
Calculated MW	114kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Sequence-specific RNA-binding protein that regulates translation and mRNA stability by binding the 3'-UTR of mRNA targets. Its interactions and tissue specificity suggest that it may be required to support proliferation and self-renewal of stem cells by regulating the translation of key transcripts.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.