

SAM68 Conjugated Antibody

Catalog No: #C56691



Package Size: #C56691-AF350 100ul #C56691-AF405 100ul #C56691-AF488 100ul #C56691-AF555 100ul #C56691-AF594 100ul #C56691-AF647 100ul #C56691-AF680 100ul #C56691-AF750 100ul #C56691-Biotin 100ul #C56691-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SAM68 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	SAM68 Antibody detects endogenous levels of total SAM68
Immunogen Description	A synthesized peptide derived from human SAM68
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Khdrbs1; p62; p68; Sam68;
Accession No.	Uniprot:Q07666
Calculated MW	68kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Recruited and tyrosine phosphorylated by several receptor systems, for example the T-cell, leptin and insulin receptors. Once phosphorylated, functions as an adapter protein in signal transduction cascades by binding to SH2 and SH3 domain-containing proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.