

Engrailed 1 Conjugated Antibody

Catalog No: #C56692

Package Size: #C56692-AF350 100ul #C56692-AF405 100ul #C56692-AF488 100ul #C56692-AF555 100ul #C56692-AF594 100ul #C56692-AF647 100ul #C56692-AF680 100ul #C56692-AF750 100ul #C56692-Biotin 100ul #C56692-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Engrailed 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Engrailed 1 Antibody detects endogenous levels of total Engrailed 1
Immunogen Description	A synthesized peptide derived from human Engrailed 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	EN1; Engrailed 1; Engrailed; Hu-En-1; Mo-en.1;
Accession No.	Uniprot:Q05925
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Belongs to the engrailed homeobox family.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.