

RAD23A Conjugated Antibody

Catalog No: #C56693

Package Size: #C56693-AF350 100ul #C56693-AF405 100ul #C56693-AF488 100ul #C56693-AF555 100ul #C56693-AF594 100ul #C56693-AF647 100ul #C56693-AF680 100ul #C56693-AF750 100ul #C56693-Biotin 100ul #C56693-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	RAD23A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	RAD23A Antibody detects endogenous levels of total RAD23A
Immunogen Description	A synthesized peptide derived from human RAD23A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	hHR23A; HR23A; RAD 23a; RAD23A;
Accession No.	Uniprot:P54725
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Multiubiquitin chain receptor involved in modulation of proteasomal degradation. Binds to 'Lys-48'-linked polyubiquitin chains in a length-dependent manner and with a lower affinity to 'Lys-63'-linked polyubiquitin chains.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.