

# ApoA4 Conjugated Antibody

Catalog No: #C56694

Package Size: #C56694-AF350 100ul #C56694-AF405 100ul #C56694-AF488 100ul #C56694-AF555 100ul #C56694-AF594 100ul #C56694-AF647 100ul #C56694-AF680 100ul #C56694-AF750 100ul #C56694-Biotin 100ul #C56694-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ApoA4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	ApoA4 Antibody detects endogenous levels of total ApoA4
Immunogen Description	A synthesized peptide derived from human ApoA4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Apo AIV; APOA 4; ApoA IV; APOA4; Apolipoprotein A IV; Apolipoprotein A4;
Accession No.	Uniprot:P06727
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

APOA4 (apolipoprotein A-IV) is a component of HDL and chylomicrons. Its primary site of synthesis is the intestine, in association with lymph chylomicron particles.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.