

UNG Conjugated Antibody

Catalog No: #C56696

Package Size: #C56696-AF350 100ul #C56696-AF405 100ul #C56696-AF488 100ul #C56696-AF555 100ul #C56696-AF594 100ul #C56696-AF647 100ul #C56696-AF680 100ul #C56696-AF750 100ul #C56696-Biotin 100ul #C56696-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	UNG Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	UNG Antibody detects endogenous levels of total UNG
Immunogen Description	A synthesized peptide derived from human UNG
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	UNG; DGU; HIGM4; HIGM5; Uracil-DNA glycosylase; UNG1; UNG2; UDG; UNG15;
Accession No.	Uniprot:P13051
Calculated MW	32kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Excises uracil residues from the DNA which can arise as a result of misincorporation of dUMP residues by DNA polymerase or due to deamination of cytosine.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.