

# ADA Conjugated Antibody

Catalog No: #C56699

Package Size: #C56699-AF350 100ul #C56699-AF405 100ul #C56699-AF488 100ul #C56699-AF555 100ul #C56699-AF594 100ul #C56699-AF647 100ul #C56699-AF680 100ul #C56699-AF750 100ul #C56699-Biotin 100ul #C56699-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ADA Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Rat
Specificity	ADA Antibody detects endogenous levels of total ADA
Immunogen Description	A synthesized peptide derived from human ADA
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ada; ADA1;
Accession No.	Uniprot:P00813
Calculated MW	41kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Catalyzes the hydrolytic deamination of adenosine and 2-deoxyadenosine. Plays an important role in purine metabolism and in adenosine homeostasis. Modulates signaling by extracellular adenosine, and so contributes indirectly to cellular signaling events.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.