

# MAGE 1 Conjugated Antibody

Catalog No: #C56701

Package Size: #C56701-AF350 100ul #C56701-AF405 100ul #C56701-AF488 100ul #C56701-AF555 100ul #C56701-AF594 100ul #C56701-AF647 100ul #C56701-AF680 100ul #C56701-AF750 100ul #C56701-Biotin 100ul #C56701-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MAGE 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	MAGE 1 Antibody detects endogenous levels of total MAGE 1
Immunogen Description	A synthesized peptide derived from human MAGE 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CT1.1; MAGE1; MAGE1A; MAGEA1;
Accession No.	Uniprot:P43355
Calculated MW	46kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

May be involved in transcriptional regulation through interaction with SNW1 and recruiting histone deacetylase HDAC1. May inhibit notch intracellular domain (NICD) transactivation. May play a role in embryonal development and tumor transformation or aspects of tumor progression.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.