

HLA B7 Conjugated Antibody

Catalog No: #C56703

Package Size: #C56703-AF350 100ul #C56703-AF405 100ul #C56703-AF488 100ul #C56703-AF555 100ul #C56703-AF594 100ul #C56703-AF647 100ul #C56703-AF680 100ul #C56703-AF750 100ul #C56703-Biotin 100ul #C56703-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	HLA B7 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	HLA B7 Antibody detects endogenous levels of total HLA B7
Immunogen Description	A synthesized peptide derived from human HLA B7
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HLA B; HLAB;
Accession No.	Uniprot:P01889
Calculated MW	40kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

HLA B belongs to the HLA class I heavy chain paralogues. This class I molecule is a heterodimer consisting of a heavy chain and a light chain (beta 2 microglobulin). The heavy chain is anchored in the membrane.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.