

MYF5 Conjugated Antibody

Catalog No: #C56708

Package Size: #C56708-AF350 100ul #C56708-AF405 100ul #C56708-AF488 100ul #C56708-AF555 100ul #C56708-AF594 100ul #C56708-AF647 100ul #C56708-AF680 100ul #C56708-AF750 100ul #C56708-Biotin 100ul #C56708-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MYF5 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	MYF5 Antibody detects endogenous levels of total MYF5
Immunogen Description	A synthesized peptide derived from human MYF5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Myf5; Class C basic helix loop helix protein 2
Accession No.	Uniprot:P13349
Calculated MW	37kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.