

GGT1 Conjugated Antibody

Catalog No: #C56711

Package Size: #C56711-AF350 100ul #C56711-AF405 100ul #C56711-AF488 100ul #C56711-AF555 100ul #C56711-AF594 100ul #C56711-AF647 100ul #C56711-AF680 100ul #C56711-AF750 100ul #C56711-Biotin 100ul #C56711-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GGT1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	GGT1 Antibody detects endogenous levels of total GGT1
Immunogen Description	A synthesized peptide derived from human GGT1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GGT; GTG; CD224; GGT 1;
Accession No.	Uniprot:P19440
Calculated MW	50-60kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Initiates extracellular glutathione (GSH) breakdown, provides cells with a local cysteine supply and contributes to maintain intracellular GSH level. It is part of the cell antioxidant defense mechanism. Catalyzes the transfer of the glutamyl moiety of glutathione to amino acids and dipeptide acceptors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.