

SLAMF1 Conjugated Antibody

Catalog No: #C56713

Package Size: #C56713-AF350 100ul #C56713-AF405 100ul #C56713-AF488 100ul #C56713-AF555 100ul #C56713-AF594 100ul #C56713-AF647 100ul #C56713-AF680 100ul #C56713-AF750 100ul #C56713-Biotin 100ul #C56713-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	SLAMF1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	SLAMF1 Antibody detects endogenous levels of total SLAMF1
Immunogen Description	A synthesized peptide derived from human SLAMF1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD150; CDw150; SLAM; SLAMF1;
Accession No.	Uniprot:Q13291
Calculated MW	85kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

High-affinity self-ligand important in bidirectional T-cell to B-cell stimulation. SLAM-induced signal-transduction events in T-lymphocytes are different from those in B-cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.