

# RCC1 Conjugated Antibody

Catalog No: #C56714

Package Size: #C56714-AF350 100ul #C56714-AF405 100ul #C56714-AF488 100ul #C56714-AF555 100ul #C56714-AF594 100ul #C56714-AF647 100ul #C56714-AF680 100ul #C56714-AF750 100ul #C56714-Biotin 100ul #C56714-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RCC1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	RCC1 Antibody detects endogenous levels of total RCC1
Immunogen Description	A synthesized peptide derived from human RCC1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	RCC1; RCC1-I; SNHG3-RCC1; CHC1;
Accession No.	Uniprot:P18754
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Guanine-nucleotide releasing factor that promotes the exchange of Ran-bound GDP by GTP. Involved in the regulation of onset of chromosome condensation in the S phase. Binds both to the nucleosomes and double-stranded DNA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.